

**CONTINUUM OF CARE
FOR TUBE FEEDING
KIDDOS AND ADULTS**

INCUBATOR

The FreeArm with Peg end holds a variety of syringe sizes for feeding and venting inside the GE Giraffe OmniBed.

CLINICAL VIDEO



CHAIR

Change to the clamp end for the FreeArm to hold the syringe during kangaroo care. No need to tape or rig the syringe. The FreeArm clamps to flat surfaces and rounded bars.

OPEN CRIB OR BED

With the clamp end of the FreeArm, the FreeArm continues to hold syringes for feeding and venting in an open crib or on a hospital bed, at the same height/rate, for consistent feeds.





SYRINGE FEEDS

The FreeArm can be discharged home with the tube feeding family to continue feeding at the same height/rate as in the hospital. DMES also provide the FreeArm and bill insurance.

[HOME VIDEO](#)

PUMP FEEDS

The FreeArm supports lightweight pumps up to 2.3 lbs and feedbags up to 1,200mL. No need to drag an IV pole upstairs, downstairs, in the car, on an airplane, etc. Making it easier to maintain nutritional compliance.

[ABSTRACT](#)



EASIER MEALS

The FreeArm works well with other products that make tube feeding easier and less messy. The Bolee Bag is a great example. Easier eating on-the-go improves quality of life for tube feeding families.

[BOLEE BAG](#)



FREEARMCARE.COM

INFO@FREEARMCARE.COM